

TAILORED FOR YOUR LIFESTYLE

This is your life and your journey. Choose how you want it to be.



CHOOSE YOUR DRIVE WHEEL STYLE

Flexibility is key - three different wheel styles open up your world like never before. Easily interchangeable to suit your daily terrain.



WHEEL STYLE	URBAN 14"	CROSSOVER 14"	OFF-ROAD 14"
WIDTH OF BASE	610mm / 24"	635mm / 25"	660mm / 26"
TERRAIN	For tight indoors and life around town	For everything in between	For going further off-road

CHOOSE YOUR COLOUR STYLE


Your chair should reflect your style and personality. Choose from a solid metallic accent colour or a hydro-dipped pattern.



MAGIC 360 SPECIFICATIONS

USER WEIGHT LIMIT	160kg / 350lb	
MAX. SPEED	10kph / 6mph	
MAX. SPEED (POWER LIFT)	5kph / 3mph	
LENGTH OF BASE	983mm / 38.7"	
WIDTH OF BASE	610mm / 24"	635mm / 25"
SEAT TO FLOOR HEIGHT	435mm / 17"	
GROUND CLEARANCE	90mm / 3.5"	
CASTORS (front & back)	Ø200mm / 8" solid	
BASE TURNING RADIUS	555mm / 21"	
MOTORS	2 x 800W gear in line	
CHAIR WEIGHT (minimum)	base, batteries, rehab seating, centre post legrest, no power lift functions	149kg / 328lb
(standard)	base, batteries, rehab seating with power lift, tilt, recline and legrest	175kg / 386lb
DRIVE ELECTRONICS	R-net	
BATTERY SIZE	73Ah gel	
SEATING (depth & width)	305-560mm / 12-22"	

STANDARD POWER FUNCTIONS



Magic 360
GO YOUR OWN WAY



Please remember that your wheelchair has electric motors and must NOT be driven through water, into rivers, creeks or the sea. Check online for current version. The information and specifications contained in this brochure are correct and up-to-date at the time of printing and may be subject to change. ©Magic Mobility. Product of Australia.

Magic Mobility +61 3 8791 5600 Email enquiries@magicmobility.com.au

Explore with us www.magicmobility.com.au



Scan QR for agent locations



M360_ENG/UY24_V07





NO JIGGLE

Improved all-wheel suspension
Engineered by the all terrain experts, and with fully independent all-wheel suspension to take care of you and your body.



KEEP UP THE PACE

5kph/3mph with up to 12" lift
Keep moving even when elevated, because the world doesn't stop for anyone.

GO ANYWHERE VERSATILITY

Choose from three interchangeable wheel options
Go where you choose, from tight indoors and offices to parks, gardens and off-road terrains.



Urban 14" drive wheels
for ultimate manoeuvrability (610mm / 24" wide)



Crossover 14" drive wheels
for added suspension and comfort (635mm / 25" wide)



Off-road 14" drive wheels
for ultimate traction outdoors (660mm / 26" wide)

CONFIGURE YOUR CHAIR TO SUIT WHERE YOU NEED TO GO AND

WHAT YOU WANT TO DO

360° Freedom

The Magic 360 is the world's most versatile compact chair with off-road capacity. Tailored for the world you live in, the Magic 360 lets you choose the life you want. Configure the body style of the chair to suit your tastes. Select the wheels for the adventures you choose. Tailor the world to suit you. Your chair. Your terrain. Your choice. Because the world's not flat™

FROM TIGHT INDOORS, URBAN STREETSCAPES TO OFF-ROAD

URBAN MANOEUVRABILITY NOW GOES OFF-ROAD

610-660mm wide compact body
The first chair that's compact enough to move easily indoors, and equipped with everything you need for the great outdoors.



2X THE POWER TO PERFORM

800 watt motors
Only the most powerful motors will enable you to take on more in life with absolute confidence.