

Name of Standards	Standard No.	Extreme X8	Frontier V6	Frontier V4 RWD	Frontier V4 FWD	Magic 360	XT2	XT4
Wheelchairs - Part 1: Determination of Static Stability	ISO 7176-1	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 2: Determination of Dynamic Stability of electrically powered wheelchairs	ISO 7176-2	✓	✓	AS /NZS 3696.2	AS /NZS 3696.2	✓	✓	✓
Wheelchairs - Part 3: Determination of the effectiveness of brakes	ISO 7176-3	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 4: Energy Consumption of electric wheelchairs and scooters for determination of theoretical distance range	ISO 7176-4	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 5: Determination of overall dimensions, mass and manoeuvring space	ISO 7176-5	✓	✓	AS /NZS 3696.5	AS /NZS 3696.5	✓	✓	✓
Wheelchairs - Part 6: Determination of maximum speed, acceleration and deceleration of electric wheelchairs	ISO 7176-6	✓	✓	AS /NZS3 696.6	AS /NZS3 696.6	✓	✓	✓
Wheelchairs - Part 7: Measurements of seating and wheel dimensions	ISO 7176-7	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 8: Requirements and test methods for static, impact and fatigue strengths	ISO 7176-8	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 9: Climatic tests for electric wheelchairs	ISO 7176-9	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 10: Determination of obstacle-climbing ability of electrically powered wheelchairs	ISO 7176-10	✓	✓	✓	✓	✓	✓	✓

Name of Standards	Standard No.	Extreme X8	Frontier V6	Frontier V4 RWD	Frontier V4 FWD	Magic 360	XT2	XT4
Wheelchairs - Part 14: Power and control systems for electrically powered wheelchairs and scooters - Requirements and test methods	ISO 7176-14	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 15: Requirements for information disclosure, documentation, and labelling	ISO 7176-15	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 16: Resistance to ignition of postural support devices	ISO 7176-16	✓	✓	✓	✓	✓	✓	✓
Wheelchairs Part 19: Wheeled mobility devices for use as seats in motor vehicles	ISO 7176-19	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 21: Requirements and test methods for electromagnetic compatibility of electrically powered wheelchairs and scooters, and battery chargers	ISO 7176-21	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 25: Batteries and chargers for powered wheelchairs	ISO 7176-25	✓	✓	✓	✓	✓	✓	✓
Wheelchairs - Part 30: Wheelchairs for changing occupant posture	ISO 7176-30	✓	✓	x	x	✓	✓	✓
Electrically powered wheelchairs, scooters, and their chargers - Requirements and test methods	EN12184	✓	✓	AS /NZS 3695.2	AS /NZS 3695.2	✓	✓	✓
Quality Management	ISO 13485 ISO 9001	✓	✓	✓	✓	✓	✓	✓

**Contact Magic Mobility if you require further information:**

3 International Court, Scoresby Vic 3179  
 Tel +61 3 8791 5600  
 Email [sales@magicmobility.com.au](mailto:sales@magicmobility.com.au)  
<http://www.magicmobility.com.au>